

Towards Riemannian EEG classifiers to detect awake and anesthetized states using Median Nerve Stimulation

Valérie Marissens Cueva

valerie.marissens@loria.fr

Supervisors:

Fabien Lotte, Laurent Bougrain, Sébastien Rimbart

Accidental Awareness during General Anesthesia (AAGA)




Yuko Shimizu

Monitors to prevent AAGA

- **Awareness still occurs** with current monitors [6, 7, 8, 9, 10].

Reasons:

- Indices were developed from adult cohorts [11]
- Indices might be falsely increased with EMG activity [7, 8] and are less accurate with age [12]



Bispectral Index monitor, Sigl and Chamoun, J Clin Monit, 1994

[6] Myles et al. Lancet, 2004.

[7] Sebel et al. Anesth. Analg. 2004.

[8] Bruhn et al. Anesth, 2000.

[9] Aho et al. Acta Anaesthesiol Scan, 2012.

[10] Schneider et al. Br J Anaesth, 2003.

[11] Bruhn et al. Br J Anaesth, 2006.

[12] Ni et al. Br J Anaesth, 2019.

Long-term objective

Propose a **new Brain-Computer Interface (BCI)** based on motor imagery to:

- monitor the motor cortex with EEG
- **detect movement intentions (MI)**
- and alert the medical staff [13]



feasible with current BCIs, although the MI vs rest classification remains low [14] for a reliable device to be used in a hospital

[13] [Rimbert](#) et al, Front. Neurosci, 2019.

[14] Lotte et al, J. Neural Eng, 2018

Median Nerve Stimulation (MNS)-based BCI

- **Painless stimulation** of the median nerve, with an intensity between 8 and 15 mA.



Illustration of the regular branching of the median nerve [15]






Illustration of ERD/ERS induced by MI, MNS and MNS during MI, in the alpha and beta frequency bands [13].

[13] Rimbart et al, Front. Neurosci, 2019.

[15] Lanz, J. of Hand Surgery, 1977.

Median Nerve Stimulation (MNS)-based BCI

- **Painless stimulation** of the median nerve, with an intensity between 8 and 15 mA.
- Combining motor intention (MI) and MNS significantly impacts MNS-generated patterns, leading to **better MI detection accuracy (+12%)** in patients under light sedation [13].
- Provides a trigger.



Distribution of average classification accuracies (n = 16), obtained by a Riemannian classifier
***p-value < 0.001 [13].

[13] [Rimbert et al](#), Front. Neurosci, 2019.

[15] [Lanz](#), J. of Hand Surgery, 1977.

Problematic

General anesthesia $\geq 4 \mu\text{g/ml}$ [19]

- The MNS pattern is still visible under light sedation ($\leq 1 \mu\text{g/ml}$) [16, 17].
- We are still able to detect it at deeper concentrations of anesthesia [18].

Is it possible to use median nerve stimulation
to detect recovery of consciousness
under general anesthesia?

[16] [Guerci](#) et al, SFAR, 2021.

[17] [Rimbert](#) et al, IEEE SMC, 2019.

[18] [Rimbert](#) et al, ASA Annual Meeting, 2023.


[19] Anderson et al, Anesthesiology, 2019.

Clinical protocol at CHU Brugmann

12 patients: 6 females; 50 ± 14.3 years old




[20] [Rimbert](#) et al. JMIR Res. Protoc, 2023.



Clinical protocol at CHU Brugmann


Surgeries:

- Anal fistula (40 min)
- Cholecystectomy (40 min)
- Left laparoscopic salpingo-oophorectomy (1h15)
- Anal sphincteroplasty (1h20)
- Laparoscopy for endometriosis (2h10)
- Transurethral resection of the prostate (1h40)
- Hemorrhoids (40 min)
- Ileostomy closure (1h30)
- Orthopedic plate change (2h20)



Data collection: protocol

MNS: Median Nerve Stimulation




- 1 run = 150 stimulations spaced by 3 to 4 seconds
- the duration and number of runs depended on the length of the surgery

Can we differentiate
awake vs anesthetized
states?

Impact of GA on the MNS patterns

PSR: Post-Stimulation Rebound
ERD: Event Related Desynchronization
PMBR: Post-Movement Beta Rebound



Grand average time frequency analysis across 8 subjects for electrode C4

MNS sensorimotor modulations [21, 22] are strongly **influenced by general anesthesia**: PSR, ERD and PMBR disappear


[21] [Rimbert](#) et al. 2th IEEE Intern. Winter Conf. on BCI, 2024.

[22] [Marissens Cueva](#) et al. J. Neural Eng, *in submission*, 2024.


Riemannian Minimum Distance to the Mean (MDM)



MNS-awake




MNS-anesthesia



[23] Congedo et al. 8th Inter. BCI Conf. 2019.

[24] Barachant et al. IEEE Trans. Biomed. Eng. 2012.


Riemannian Minimum Distance to the Mean (MDM)



[23] Congedo et al. 8th Inter. BCI Conf. 2019.

[24] Barachant et al. IEEE Trans. Biomed. Eng. 2012.


Riemannian Minimum Distance to the Mean (MDM)




[23] Congedo et al. 8th Inter. BCI Conf. 2019.

[24] Barachant et al. IEEE Trans. Biomed. Eng. 2012.

Offline classification: MNS-awake vs MNS-anesthesia



Test set balanced accuracy of a Riemannian Minimum Distance to the Mean [23, 24] (10 subjects). The orange line indicates the average test set balanced accuracy (97%).





Problem Classifier trained with **deep anesthesia trials**, which are **not available before the surgery** to calibrate the BCI

[23] Congedo et al. 8th Inter. BCI Conf. 2019.

[24] Barachant et al. IEEE Trans. Biomed. Eng. 2012.



One-Class Riemannian Minimum Distance to the Mean (MDM)



New covariance matrix P of class k :

$K = C_1$ if $\delta_R(P, P_i) \leq \delta_{th}$ else C_2
with δ_R the Riemannian distance:
 $\delta_R(P, P_i) = \|\log(P^{-1/2} P_i P^{-1/2})\|_F$


One-Class Riemannian Minimum Distance to the Mean (MDM)



New covariance matrix $\mathbf{P}$ of class k :



$\mathbf{K} = \mathbf{C}_1$ if $\delta_R(\mathbf{P}, \mathbf{P}_i) \leq \delta_{th}$ else $\mathbf{C}_2$
with δ_R the Riemannian distance:
 $\delta_R(\mathbf{P}, \mathbf{P}_i) = \|\log(\mathbf{P}^{-1/2} \mathbf{P}_i \mathbf{P}^{-1/2})\|_F$

One-Class classification trained on MNS-awake




Test set balanced accuracy of a **One-Class Riemannian Minimum Distance to the Mean** [25] (12 subjects). The orange line indicates the average test set balanced accuracy (85%).

Distance between covariance matrices and MNS-awake class centroid (C)



Subject 8 (97.3 %)


Distance between covariance matrices and MNS-awake class centroid (C_j)



Conclusions

- It is indeed possible to use **median nerve stimulation** to detect the recovery of consciousness under general anesthesia
- A One-Class Riemannian Minimum Distance to the Mean algorithm **trained with only the awake data** correctly differentiates between MNS-awake and MNS-anesthesia (test balanced accuracy of **85.44 %**)

Perspectives for a new MNS-based BCI to detect AAGA



Acknowledgments



Fabien Lotte
Research Director
Inria Center at the University of Bordeaux



Laurent Bougrain
Associate Professor
Université de Lorraine



Sébastien Rimbart
Inria Starting Faculty Position
Inria Center at the University of Bordeaux



Seyed Javad Bidgoli



Denis Schmartz



Claude Meistelman



Ana Maria Cebolla
Alvarez



Guy Cheron



Philippe Guerci



Thank you for your attention!

Contact: valerie.marissens@loria.fr



Monitors to prevent AAGA



NeuroSENSE monitor

www.neurowavesystems.com




SedLine monitor

www.maasimo.fr



Narcotrend monitor

www.narcotrend.de



Bispectral Index monitor, Sigl and Chamoun, J Clin Monit, 1994


www.nursekey.com

Threshold selection

Δ : distances between training MNS-awake set covariance matrices and MNS-awake centroid


1. $\text{mean}_{\Delta} + 3\text{std}_{\Delta}$
2. $\text{median}_{\Delta} + 3\text{std}_{\Delta}$
3. $\text{median}_{\Delta} + 3\text{mad}_{\Delta}$

with std: standard deviation and mad: median absolute deviation




Which training class?

- i. 1st preoperative run
- ii. 1st intraoperative run
- iii. 2nd intraoperative run




Electrocautery-related EEG artifacts



OCMDM


- 250:1000 ms
- mu+beta frequency band
- all 64 electrodes




	condition i	condition ii	condition iii
Train	65:296	108:150	70:199
Test intraop	414:1480	306:1288	267:1366
Test preop	28:150	94:445	94:445

Number of trials used for the 13 subjects (9 for condition ii)
(min:max)


Classification results per subject




Classification results per subject



Classification results per subject



Informations des patients



Offline classification: MNS-awake vs MNS-anesthesia

- Classifier: Riemannian Minimum Distance to the Mean
- 64 electrodes
- Preprocessing :
 - downsampling 128 Hz
 - 8-30 Hz band pass filter
 - epochs : 250 to 1000 ms